


# EXHAUST ELEMENTS

COMPLEMENTS

The background features a central vertical bar in a light grey color. This bar is flanked by several large, overlapping triangles in shades of purple and grey. The overall composition is abstract and geometric, with sharp lines and a limited color palette.



### DIMENSIONS

W : Width  
600/900/1200/1500  
mm.

D : Depth  
550 mm.

H: Height  
600 mm.

Ventilation connection  
Ø 200 mm.

Recommended flow  
extraction  
300/400/500/600  
m<sup>3</sup>/h.

### \* TECHNICAL CHARACTERISTICS

Trapezoidal aspiration hood without the enclosure to evacuate thermal loads or not harmful steam in an open way in the laboratory.

### \* MATERIALS

Structure with possibility of manufacturing in PP or AISI 304 stainless steel.

